

DNA FACTSHEET

CONCENTRATION BY ABSORBANCE

Spectrophotometry can be used to estimate DNA or RNA concentration and analyse the purity of the preparation.

For a 1cm pathlength:

- For double stranded DNA $1A_{260} = 50\mu\text{g/mL}$
- For single stranded DNA $1A_{260} = 33\mu\text{g/mL}$

Based on extinction coefficients of DNA in water

The OD does not give any indication of the size of the DNA

PURITY BY ABSORBANCE

- Pure DNA has an A_{260}/A_{280} ratio ≥ 1.8
- An A_{260}/A_{280} ratio < 1.8 indicates contamination with proteins or aromatic substances
- An A_{260}/A_{280} ratio > 2.0 indicates a possible contamination with RNA

Note a high A_{230} reading can also indicate contaminants in the sample

MOBILITY OF DIFFERENT SIZE DNA FRAGMENTS ON AGAROSE GELS

Agarose %	Effective resolution range of linear ds DNA fragments (kb)
0.5	30 to 1
0.7	12 to 0.8
1.0	10 to 0.5
1.2	7 to 0.4
1.5	3 to 0.2

BIOLOGICAL BUFFERS

The control

CITRATE BUFFER: USEFUL RANGE PH3.0 TO 6.2

Table 1 : Citrate buffer components

Component	Name	Synonyms	Formula	Mol. Wt.
A	Citric acid monohydrate	N/A	$C_6H_8O_7 \cdot H_2O$	210.14
	Or Citric acid	N/A	$C_6H_8O_7$	192.12
B	Trisodium citrate dihydrate	Sodium citrate tribasic dihydrate, Citric acid trisodium salt dihydrate	$C_6H_5O_7Na_3 \cdot 2H_2O$	294.12

OTHER COMMON BUFFER RECIPES

REFERENCES
